BUBBLE FOOT MASSAGER
INSTRUCTION MANUAL

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Thank you for buying
an ORIGINAL Product of

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INSTRUCTION MANUAL

FEATURES:

1. Automatic heat insulation: There's a high-tech PTC heating device at the bottom of the foot spa in order to increase the heat efficiently to a certain temperature automatically and can also keep the water temperature. In this way you can enjoy the fun of foot health comfortable temperatures at well. It can be heated to a comfortable temperature from 40°C to 45°C.
2. Vibration massage: There is hundreds of massage columns at the bottom of the foot spa, put your feet on them after the vibration so you can stimulate the acupuncture points completely
3. Air Bubble Impact: There are strong air bubbles from 2 surfs around both sides of the bottom which could massage you and reduce fatigue, soften skin to reach the effect of therapeutic massage.
4. Infrared hot moxibustion: The infrared device can be ad hoc foot full of hot moxibustion, it can activate collaterals, and promote blood circulation.
5. Aromatherapy Conditioner: you can put all kinds of drops of essential oils in massage basin and at the same time you may have the full enjoyment and refreshing quickly relieve fatigue.
6. Foldable body design: With foldable body to easy stock. You can fold the body with especial material.

IMPORTANT SAFETY INSTRUCTIONS:

When using an electrical appliance, basic precautions should always be followed, including the following:

READ ALL INSTRUCTIONS BEFORE USING

DANGER – to reduce the risk of electric shock

- Always unplug the appliance from the electrical outlet immediately after using and before cleaning.
- DO NOT reach for an appliance that has fallen into water. Unplug it immediately.
- DO NOT use while bathing or showering.
- DO NOT place or store appliance where it can fall or be pulled into a tub or sink. DO NOT place in or drop into water or other liquid.
- DO NOT attempt to plug or unplug unit while feet in the water.

WARNING:

To reduce the risk of burns, fire, electric shock or injury to persons:

- This appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the appliance by a person responsible for their safety.
- Children should be supervised to ensure that they do not play with the appliance.
- An appliance should never be left unattended when plugged in. Unplug from outlet when not in use before putting on or taking off parts or attachments.
- Close supervision is necessary when this appliance is used by, on, or near children, invalids or disable persons.
- Use this appliance only for its intended use as described in this manual. Do not use attachments not recommended, specifically any attachments not provided with this unit.
- Never operate this appliance with a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or dropped into water. Return it to the service centre for examination and repair.
- Keep cord away from hot surfaces.
- Never block the openings of the appliance, or place it on a soft surface, such as a bed or a couch, where the air openings may be blocked. Keep air openings free of lint, hair, dust etc.
- Never use while sleeping or fall asleep while using the massager.
- Never drop or insert any object into any opening.
- Do not operate where aerosol (spray) products are being used or where oxygen is being administered.
- Do not operated under blanket or pillow. Excessive heating can occur and cause fire, electrical shock or injury to persons.
- Do not carry this appliance by the supply cord or use the cord as a handle.
- To disconnect, turn all controls to position "OFF", then remove the plug from outlet.
- Do not attempt to stand on or in foot massager. Use only while seated.
- Use heated surface carefully. May cause serious burns. Do not use over sensitive skin areas or in the presence of poor circulation. The unattended use of heat by children or incapacitated persons may be dangerous.
- Not use outdoors.

- Connect this appliance to a properly grounded outlet only. See grounding instructions.

Caution:

All serving of this foot massager must be performed by authorized service personal only.

SAVE THIS INSTRUCTIONS! FOR HOUSEHOLD USE ONLY!

Caution – Please read all instructions carefully before operating:

- If you have any concerns regarding your health, consult your doctor before using foot massager.
- Individuals with pacemakers should consult a physician before using this product.
- In case of pregnancy, diabetes or illness, consult your doctor before using.
- We recommended using this unit for no longer than 30min at a time. Extensive use could cause the product to overheat. Should this occur, discontinue use and allow the unit to cool before operating again.
- If you experience pain in muscles or joint for a prolonged period of time, discontinue use and consult your doctor. Persistent pain could be a symptom of a more serious condition.
- Never use any massager on open wounds, discoloured areas, or any area of the body that is swollen, burned, inflamed or where skin eruptions or sores are present.
- Massage should be pleasant and comfortable. Should pain or discomfort result, discontinue use and consult your physician.

SPECIAL WARNINGS:

- Never immerse foot massager in water or other liquids.
- Always check the temperature of the water before placing your feet into foot massager.
- Do not stand in the unit. Sit down and place your feet in foot massager.
- Always unplug the unit when not in use and when filling, emptying, cleaning or moving the unit.
- Ensure that your hands dry when operating the switch or removing the plug.
- Do not leave foot massager unattended while it is operating.

INSTRUCTIONS FOR USE:

Ensure that the unit is unplugged before filling with water. **Do not fill water over the fill line.** Add water to the basin first according to the required height of water level and turn on power, use only the 230V power supply, infrared lights, press the power button, indicator light turns red. (Non-water use)

1. The first button for the power key.
2. The second button for the vibration/ bubble/ heating function
3. To empty: first unplug the unit. Tip the unit allowing water to spill out along the drain spout.

Never pour water over the control switch.

MAINTENANCE:

TO CLEAN:

Unplug the unit and allow it to cool before cleaning. Clean only with a soft, damp cloth. Never allow water or any other liquid to come into contact with the control switch.

- After cleaning, wipe surfaces with a dry cloth.
- Never use abrasives, detergents or cleaning agents.
- Periodically straighten cord if it becomes twisted.

TO STORE:

Unplug the unit and allow it to cool.

- Store in its box or a clean, dry place.
- To avoid breakage, do not wrap the power cord around the unit.
- Do not hang the unit by the power cord.

SPECIFICATIONS:

220-240V • 50/60Hz • 450W

Correct Disposal of this product

You can help protect the environment! Please remember to respect the local regulations: hand in the non-working electrical equipments to an appropriate waste disposal center..

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